

## **ANALYZER SOLUTIONS FOR YOUR PROCESS!**

# **AMETEK ProLine Process Mass Spectrometer**

Fast and reliable multiport, multicomponent process analysis

## PROCESS MS AT A PROCESS GC PRICE

Offering high performance without the high price, the ProLine Series of process mass spectrometers from AMETEK Process Instruments provides features and a level of performance associated with much more expensive mass spectrometer systems. At last, you can afford to use MS for applications currently using GC, NIR, FTIR or multiple single-component analytical devices.

#### FIELD-PROVEN RELIABILITY

AMETEK has more than 20 years of experience in mass spectrometer design and manufacture, with many thousands of analyzers already operating reliably. The field-proven analyzer components, combined with AMETEK's more than 30 years' experience in the process analyzer business, means you can depend on ProLine to take care of your critical process control requirements.

# **AUTOMATIC OPERATION**

Automatic calibration and verification maintains the performance of the analyzer and allows it to quantify the amount of gas species in multicomponent samples. Operation of the ProLine is fully automatic, providing unattended analysis of process streams and calibration gases. External process parameters can be used to trigger events. User-definable data presentation lets you set up the displays and outputs for your application needs.

# MULTIPOINT, MULTICOMPONENT ANALYSIS

Real-time process monitoring of up to 32 components is straightforward with ProLine's powerful easy-to-use software. Complex overlapping spectra are handled automatically. The ProLine provides data output directly in concentration units (PPM, %, etc.) with high sample rate and data throughput.

Reactor 1202	
Nitrogen	76.478
Oxygen	17.109
Argon	1.859
C02	4.554
RQ	1.015

Multicomponent analysis. Up to 32 components, 8 or 16 streams.



The ProLine comes complete with multiport automatic valve switching, inlet manifold, and turbomolecular pumping package.

## **COMPACT**

Analyzer space is always at a premium. Instrumentation has become smaller, lighter, less expensive and better, and it has just happened to process mass spectrometers. The AMETEK ProLine will easily fit on a corner of your bench or workspace.

#### APPLICATIONS

- Fermentation/bioreactor off-gas analysis
- Fuel cell research and development
- Catalytic reactions/temperature-programmed desorption
- Semiconductor process exhaust gas monitoring
- Bulk gas analysis
- Gas blending
- Glove box monitoring
- Pharmaceutical solvent drying
- Gas scrubber efficiency
- Ethylene/ethylene oxide production
- Synthesis gas
- R&D laboratories

## **SPECIFICATIONS**

Range: 1-100 AMU. Optional: 1-200 AMU, 1-300 AMU.

Sample Ports: 8- or 16-way stream selection manifold

Inlet Pressure: 2 to 20 PSIG (14 to 140 kPa); optional sub-

atmospheric pressure

**Accuracy:** Better than 0.5% of measured value for argon in

Outputs: RS-485, RS-232, optional 4-20mA **Detection Range:** From 1 PPM to 100%

Power: 115/230 VAC, 50/60 Hz, 500W

PC Requirements: Pentium® equivalent, 600MHz or better,

Windows 98, ME, NT, 2000 or XP

Enclosure: General purpose benchtop

**Dimensions:** 14 3/8" x 24" x 15 1/2" (36.5 cm. x 61 cm. x

39.37 cm.)

Weight: 80 lb. (36 kg.)

# **MODULAR, EASY-TO-SERVICE**

The compact modular design of the ProLine sets a new standard in serviceability. The mass spectrometer electronics are contained in a single, easily replaced, plug-in unit. No electronic calibrations are necessary-the electronics unit will mate with any analyzer head. An inexpensive spare analyzer head and electronics unit can be kept on the shelf and easily exchanged by your own personnel on-site. There is no need for the expensive service contract typical of most process mass spectrometers. Self-diagnostics and automatic device event logging ensure GMP compliance. Modem support is available for additional factory diagnosis and troubleshooting.

#### HIGH PERFORMANCE MASS ANALYZER


At the heart of this versatile mass spectrometer is the proven AMETEK quadrupole mass filter coupled with a contaminantresistant enclosed ion source. The ProLine shows unequaled, remarkably clean mass separation, and delivers outstanding sensitivity and stability across the whole mass range.



Analyzer head and electronics unit.

### MULTIPORT SAMPLING SYSTEM

The multiport sampling system consists of 8 or 16 manifold-mounted valves and a capillary restriction in an enclosure heated to prevent condensation. The inlet arrangement is designed for fast clearing and offers the choice of continuous flush on any valve position. The sampling system can be configured for sampling gas pressures from 2 to 20 PSIG with an option available that allows sampling of subatmospheric pressure gases. Different calibration and analysis methods can be assigned to each sample port.



## WHY PROLINE FOR YOUR PROCESS

- Lowest price-per-sample point
- Small and compact size
- Multipoint backup for other process analyzers
- Multiport, multicomponent analysis
- Automatic sample stream switching
- Exceptional reliability
- Real-time analysis
- Scan for unwanted sample components
- High sensitivity
- Self-diagnostics

USA - Texas

Economical and versatile general troubleshooting and development tool

One of a family of innovative process analyzer solutions from AMETEK Process Instruments. Specifications subject to change without notice.

USA - Pittsburgh 150 Freeport Road Pittsburgh, PA 15238 Ph. 412-828-9040 Fax 412-826-0399

> CHINA Ph. 86 21 6426 7049 ax 86 21 6426 7054

USA - Delaware 455 Corporate Blvd. Newark, DE 19702 Ph. 302-456-4400 Fax 302-456-4444

> FRANCE Ph. 33 1 30 68 69 20 Fax 33 1 30 68 69 29



Ph. 49 21 59 91 36 0 Fax 49 21 59 91 3680

**GERMANY** 

MIDDLE EAST - Dubai Ph. 971-4-881 2052 Fax 971-4-881 2053

9750 Whithorn Drive

Houston, TX 77095

Ph. 281-463-2820

Fax 281-463-2701

CANADA 2876 Sunridge Way N.E. Calgary, AB T1Y 7H9 Ph. 403-235-8400 Fax 403-248-3550

SINGAPORE Ph. 65-484-2388 Fax 65-481-6588

